

## CHEMISTRY EXPERIMENT THEORETICAL YIELD

In this experiment sodium hydrogen carbonate (baking soda) will be reacted with hydrochloric acid to obtain a high yield of sodium chloride. The unbalanced equation for the reaction is:



You will:

- completely react an accurately measured quantity of NaHCO<sub>3</sub> and dilute HCl.
- isolate the NaCl from the other products and determine its mass.

The **theoretical yield** can be calculated by using the mole ratio from the balanced equation for the reaction.

After completing the lab, the **percent yield** can then be determined by comparing the experimental yield to the previously established theoretical yield.

### **PROCEDURE:**

#### **Day One:**

1. Determine the mass of a 100mL beaker to the nearest 0.01 gram.
2. Weigh out, precisely, your assigned mass (\_\_\_\_\_g) of sodium hydrogen carbonate into the beaker. Use a microspatula to adjust amount. Record the total mass of beaker plus contents.
3. Cover the beaker with a watch glass, convex side down and slightly off center so that the lip of the beaker is uncovered. Then:
  - Add dilute (3M) HCl down the lip of the beaker to the sodium hydrogen carbonate in the beaker.
  - Continue this process until no more reaction takes place when a drop of acid is added (the bubbling stops).
  - Gently swirl the beaker so that all of the solid contacts the acid. Do **NOT** add excess hydrochloric acid!
4. Carefully rinse the underside of the watch glass with a small amount of distilled water and collect the washings in the beaker. Wash the watch glass with tap water, dry, and return.
5. **Gently** heat the beaker and contents on a hot plate until the sodium chloride solution *just begins to boil*. This will drive off any excess HCl and begin the evaporation process.
6. Using beaker tongs, place the beaker in a drying oven to evaporate the water from the NaCl overnight.

#### **Day Two:**

1. Weigh the beaker and contents (now crystalline NaCl) to the nearest 0.01 gram. Record.
2. Clean, rinse, and dry the beaker. The NaCl can be washed down the drain.

## THEORETICAL YIELD Lab Write-up

**PURPOSE:** The purpose of this experiment is to obtain a high yield of sodium chloride and to compare that yield to a value obtained by doing a mass-to-mass stoichiometric calculation, thus validating this type of calculation.

### Day One CALCULATION:

Write the balanced equation for the reaction. Then calculate the **theoretical yield** of NaCl that should be obtained using your assigned mass of NaHCO<sub>3</sub>. This is a gram-to-gram **stoichiometry** calculation.

**(g of A → mol of A → mol of B → g of B)** Show factor-label method and units!!!

Equation:



Mass of 100mL beaker \_\_\_\_\_ grams

= Mass of NaHCO<sub>3</sub> \_\_\_\_\_ grams

Mass of beaker + NaHCO<sub>3</sub> \_\_\_\_\_ grams

*Calculation of Theoretical Yield:*

$$\underline{\hspace{1cm}} \text{ g NaHCO}_3 \left( \frac{\hspace{1cm}}{\hspace{1cm}} \right) \left( \frac{\hspace{1cm}}{\hspace{1cm}} \right) \left( \frac{\hspace{1cm}}{\hspace{1cm}} \right) = \underline{\hspace{1cm}} \text{ g NaCl}$$

**Day Two DATA:** Mass of beaker + NaCl \_\_\_\_\_ grams

### ANALYSIS/CALCULATIONS:

1. Use subtraction to find the Experimental Yield of NaCl.

2. Determine the percentage yield of NaCl.

$$\% \text{ Yield NaCl} = \frac{\text{Experimental Yield NaCl}}{\text{Theoretical Yield NaCl}} \times 100 =$$

### CONCLUSION/QUESTIONS:

1. State in a sentence what your % Yield means. This is asking you to explain the numerical result to calculation #2 above, not how you got it.

2. Based on your % Yield, how accurately did your stoichiometric calculation actually predict your product? Explain any discrepancies.

3. Calculate the volume, in mL, of 3.00 M HCl that would be theoretically required to react with the grams of NaHCO<sub>3</sub> that you used in the reaction. This is a solution stoichiometry problem. g A → mol A → mol B → mL B

$$\underline{\hspace{1cm}} \text{ g NaHCO}_3 \left( \frac{\hspace{1cm}}{\hspace{1cm}} \right) \left( \frac{\hspace{1cm}}{\hspace{1cm}} \right) \left( \frac{1000 \text{ mL}}{3.0 \text{ mol HCl}} \right) = \underline{\hspace{1cm}} \text{ mL HCl}$$